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SYNTELLIX CE 0197

1. Product description

The instruments of Syntellix AG are used in accident surgery and orthopaedics. Reusable instruments are not delivered in a sterile condition and must be cleaned, serviced and autoclaved prior to every surgical use. Before they are cleaned for the first time, original packing and protective caps must be removed completely.

2. Choice and application of instruments

The instrumnts may be applied in correspondence with their intended use only. The attending physician is responsible for the choice of the instruments to be used in the surgical procedure. Prior to any surgery, the user must verify that the instruments are in a sound condition and optimally attuned to each other and, where applicable, to the implants to be inserted. The products may be used by trained personnel only. The implantation method is to be chosen based on the product-specific surgical instructions and on the current state of medical science. In addition, the separate mounting and application instructions for screw measuring devices shall apply.

3. Combination notes

Syntellix AG does not assume liability for any complications that may occur due to the combination with products supplied by other manufacturers.

4. Warnings

Applying the instruments for any other purpose than their intended use may result in damage to the instrument and, thus, in danger to the patient and/or the user. Any overstraining of the instrument may weaken the material and, consequently, cause intraoperative product failure. Hence, checks must be performed prior to and following any use in order to ensure completeness and functionality.

5. Cleaning / disinfection

The instrumts from Königsee Implantate GmbH are exclusively thermostable. Machine cleaning and disinfection as described below is preferable to the manual method.

Choice of the cleaning agent

In order to avoid corrosion damage, mildly alkaline cleaning agents ought to be used. Sufficient rinsing with deionized water is always necessary. Where a highly alkaline cleaning agent is used, a neutralization step involving deionized water and a neutralizing agent (e.g. Mucapur Z) may have to be included in the cleaning procedure.

- organic, mineral and oxidizing acids
- organic solvents (e.g. acetone, ether, alcohol, benzine)
- oxidizing substances (e.g. peroxide)
- halogens (chlorine, iodine and bromine compounds)
- aromatized, halogenated hydrocarbons

Pre-treatment

Instruments must be deposited properly during surgery in order to prevent damage. Contaminated instruments must be deposited separately already after use. Visible residues have to be removed from the product two hours after use at the latest. For this purpose, cold running water or a disinfectant solution is to be applied. The disinfectant should be free of aldehydes and have been tested for effectiveness. The product must never be deposited in physiological saline solution since prolonged contact can lead to corrosion. Only a soft brush or a soft, clean cloth may be used for manual cleaning. Cavities are to be rinsed thoroughly with a water pistol. To the extent possible, the instruments must be disassembled and hinges must be opened.

Pre-cleaning in an ultrasonic bath

Cleaning in an ultrasonic bath can be done to support the cleaning procedures performed by hand and by means of machines. It is particularly suitable in case of instruments with hinges that cannot be disassembled. The instruments must be opened and immersed in the ultrasonic bath (at 35 kHz) containing deionized water and a cleaning additive (solution in water according to the manufacturer's instructions for disinfecting solution at 30°C±3°C). Depending on the degree of contamination of the instruments, the contact time can be selected within a range from 10 to 30 min. When taking instruments with small components (nuts and/or screws) out of the ultrasonic bath, it is to be made sure that the components are removed. In general, all instruments supplied by Syntellix AG are suitable for being treated in an ultrasonic bath.

Instrument components made of silicone, however, lessen the effect of the ultrasound due to the flexibility of the material.

Manual cleaning and disinfection

Fully immerse the implants in a disinfection bath (solution in water according to the manufacturer's instructions for disinfecting solution at 30°C±3°C) for 20 minutes; make sure that the implants cannot touch each other in the bath. In case of instruments with cavities, it is recommendable to move the parts around by hand to make sure that the internal space is rinsed thoroughly.

Following the contact time, the entire external surface of the instruments - especially areas that are hard to reach (cavities, threads and hinges) - is to be brushed using a soft, non-metallic brush for at least one minute. Couplings must be treated with special care during the cleaning procedure. Cannulas are to be brushed under running water using a rounded nylon brush which must be moved back and forth at least 10 times.

Movable sliding surfaces (e.g. hinges, couplings) have to be operated under running water at least 10 times. Finally, the instruments must be rinsed using cold deionized water for at least 10 seconds.

They are to be dried using a soft, lint-free cloth. Cavities must be dried by applying clean compressed air.

The manual cleaning and disinfection procedure has been validated using MediZym (1 % solution).

Machine cleaning and disinfection

The washer-disinfector to be used must fulfill the requirements of EN ISO 15883 and include a suitable program for thermal disinfection (AO ≥ 3000).

Process sequence:

1. Deposit and/or insert the disassembled instruments in the appropriate receptacle inside the washer-disinfector. The products must not come into contact with each other. Areas that remain unwashed are to be avoided.

2. Start des Programms.

Phase	Duration	Note
Pre-cleaning	1 min	Rinsing with cold water
Cleaning	5 min	Water (55°C ± 1°C) with cleaning additive
Neutralization*	2 min	Water with neutralization additive
Rinsing	1 min	Warm water (30°C±3°C)
Thermal disinfection	5 min	≥ 93°C with deionized water without any additives
Drying	50 min	115°C

3. Take the instruments out of the washer-disinfector after the end of the program.

*. If a highly alkaline cleaner is used, neutralization using fully deionized water with a neutralizer can be included as a step in the washing process.

This procedure has been validated applying the cleaning additive neodisher® MediClean.

6. Inspection, assembly and maintenance

After cleaning, the instruments are to be cooled down to room temperature and then inspected for damage, wear, contamination and malfunctioning parts. In case of visible residues, the products must be treated once again. The user is responsible for sorting out damaged and contaminated instruments.

Multipart instruments have to be assembled. After the instruments have cooled down, instrument oil that is suitable for autoclaving and has been tested for biocompatibility is to be used to service movable parts. Plastic surfaces must not be treated using instrument care products.

The final inspection comprises (among other things):

- Visual intactness (especially in case of cutting instruments as well as spanner flats);
- Completeness in case of multipart instruments;
- Operability of movable parts;
- Proper operation

Products that have been prepared in this manner are to be deposited in the appropriate place in the implantation set. This set is to be packed in a sterile barrier system in accordance with EN ISO 11607, i.e. in a designated sterilization wrap or in a reusable sterilization container. The instruments may also be packed in disposable sterilization packages.

7. Sterilization

Cleaned instruments are to be sterilized applying a validated autoclaving method in accordance with EN ISO 17665-1 only. A three times fractionated prevacuum process is recommended.

The following parameters must be adhered to:

Sterilization temperature	Sterilization time
132°C (±1°C)	4 min

Sufficient product drying (drying time of at least 30 min) is to be ensured. Sterilization time can be extended to a maximum of to 18 min, to comply with the recommendations of the Robert Koch Institute. The user is responsible for the validation of the sterilization techniques and sterilization equipment.

8. Additional information on cleaning and sterilization

Frequent treatment affects the instruments only slightly as long as the chemicals and rinsing water used are in accordance with the requirements described and as long as the instruments are sufficiently serviced. It must be made sure that cleaning additive residues cannot be diverted as a result of insufficient rinsing and/ or neutralization. Under certain circumstances, such residues may attack the surfaces of instruments or sterilization trays and lead to corrosion and discoloration. Products can be damaged due to the mechanical forces acting inside the washer-disinfector. Checking the instruments after cleaning is absolutely necessary. Negative effects on the functionality and quality of the products due to repeated sterilization have not been observed. The end of the product life of an instrument is normally determined by wear and damage due to use.

9. Storage conditions for instruments

Brand new instruments must be stored in dry conditions (humidity <60% and at room temperature (<30°C). The storage parameters for sterilized implants depend on the individual requirements set by the particular medical institution.

10. Repair and disposal

Faulty or incomplete instruments are to be returned to the manufacturer for repair. If repair work is carried through by unauthorized individuals, any and all warranty and guarantee claims will expire. Irreparable instruments have to be disposed of properly. The respective valid rules and regulations stipulated by the hospital operator shall apply.

11. Handling of return shipments

Returned goods can be accepted only when they are in a hygienically safe condition. The General Terms and Conditions of Syntellix AG shall apply.


12. Material and marking

The instruments are primarily made of stainless steel in accordance with EN ISO 7153-1 and can be combined with the following components:

- Surface-coated light alloys such as anodized aluminum and anodized titanium
- Polymers such as PPSU, PEEK, silicone, HGW and POM-C

The instruments are identified by way of laser marking and color coding. In correspondence with Directive 93/42 EEC, the instruments bear the CE mark. Powered instrumts as well as instrumts with a measuring function additionally show the identification number of the notified body.

All materials and marks used require careful handling and have been attuned to the treatment and sterilization methods described above. Special procedures are not necessary.

	Name and address of the manufacturer
	CE label
	"Non-sterile"
	Instructions for use must be observed.
	Detailed material specification
	Item number
	Lot number
	Quantity

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