


MAGNEZIX[®] CS 3.2

PRODUCT OVERVIEW INSTRUMENTS



INSTRUMENTS*

CS 3.2



Art. No.	Description
6032.108 1a	Screwdriver T8, One-Piece Handle Ø 1.3 mm cannulated, consisting of: 9099.001 One-Piece Handle for Screwdriver: 9032.015 Screwdriver blade T8
6032.208 1b	Screwdriver T8, Multi-Part Handle Ø 1.3 mm cannulated, consisting of: 9099.002 Multi-Part Handle for Screwdriver 9032.015 Screwdriver blade T8
9032.020 2	Drill Bit Ø 2.5/1.3 mm, cannulated, length 160/135 mm, for quick coupling
9032.021 3	Countersink Ø 3.5/1.3 mm, cannulated, for quick coupling
9032.033 4	Double Drill Guide, Ø 3.5/2.5 mm
9032.034 5	Drill Guide, Ø 2.5/1.3 mm
9032.040 6a	Guide Wire Ø 1.2 mm, with trocar tip, length 150 mm (do not reuse)
9032.041 6b	Guide Wire Ø 1.2 mm, with threaded tip, length 150 mm (do not reuse)
9032.042 7	Measuring Device, for Guide Wires Ø 1.2 mm, Guide Wire length 150 mm
9032.050 8	Cleaning Stylet Ø 1.25 mm, for Ø 1.3 mm cannulated instruments
Not shown:	
8032.001	Sterilizing Tray for MAGNEZIX® CS Ø 3.2 mm, without contents
8032.002	Lid for Sterilizing Tray, for MAGNEZIX® CS Ø 3.2 mm

*The figures are not to scale.